

# Bode SecurSwab 2

A DNA evidence collection swab system designed for the collection, protection, preservation, transport, and storage of biological samples.



FS921 and FS922

## SPECIFICATIONS

### Construction

- Integrated swab and transport device design
- Built-in dual desiccant drying system
- Available with or without barcode labeling

### Performance

- Collects biological and DNA evidence
- Integrated drying system helps preserve sample integrity
- Designed to reduce contamination and sample loss
- Reported to provide increased recoverable DNA yield
- Suitable for storage and transport of wet samples

### Quantity

- 100 per case

## Variations

Model No.	Description
FS921	Bode SecurSwab 2
FS922	Bode SecurSwab 2 with Barcodes

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI