

Swab Boxes

These 20 pt. chipboard swab boxes are designed for the collection and storage of DNA evidence such as blood and other biological fluids. They allow continued air drying after packaging and are available in multiple sizes for different swab types.



FS925



FS926



FS927

SPECIFICATIONS

Construction

- 20 pt. thick chipboard material
- .25" ventilation holes for airflow
- Available in standard, small, and suspended designs

Performance

- Stores and protects DNA swab samples
- Allows airflow for continued drying after collection
- Prevents contamination and sample degradation
- Accommodates different swab sizes for flexibility

Materials

- Chipboard

Quantity

- 100 per package

Variations

Model No.	Dimensions
FS925	3" x .5" x .5"
FS926	6" x .5" x .5"
FS927	7" x 1.5" x .875"

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI