

# Swabs- Cotton-Tipped Applicators, 3” or 6”

These cotton-tipped applicators are designed for the collection of wet and dry body fluid samples from a variety of surfaces. Available in 3” and 6” lengths, they offer versatility for different evidence collection situations.



## SPECIFICATIONS

### Construction

- Cotton-tipped applicator ends
- Wooden shaft
- Available in 3” and 6” lengths

### Performance

- Collects wet and dry body fluid samples
- Can be moistened to lift dried evidence
- Provides control or extended reach depending on length
- Suitable for a variety of surface types

### Materials

- Cotton
- Wood

### Quantity

- 2 swabs per pack, 100 packs per case

## Variations

Model No.	Length
FS929	3”
FS930	6”

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI