

Trace Evidence Collection Forensic-Swab

These forensic-swab trace evidence collection swabs are designed for collecting DNA samples and trace evidence while minimizing contamination risks. They feature a self-drying transport tube and are sterilized to ensure DNA-free integrity.



FS931

SPECIFICATIONS

Construction

- Swab with sealed transport tube
- Integrated ventilation membrane for drying
- Individually wrapped

Performance

- Collects DNA and trace evidence samples
- Self-dries within tube to prevent contamination
- Preserves sample quality for accurate analysis

Quantity

- 25 per pack

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI