

# Sterile Omni-Swabs

These sterile omni-swabs are non-invasive collection devices designed for obtaining saliva and buccal cell samples. They feature a removable swab head that allows easy transfer of collected material into tubes or plates.



FS932

## SPECIFICATIONS

### Construction

- Brush-like swab head made of absorbent filter paper
- Ejectable swab head design for sample transfer
- Non-reactive materials
- Individually sterile wrapped packaging

### Performance

- Collects saliva and buccal cell samples non-invasively
- Enables easy transfer of samples into tubes or multi-well plates
- Supports clean and efficient sample handling

### Quantity

- 100 per box

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI