

Fab-Swab Sterile DNA Controlled Standard Cotton Tipped Applicator

The fab-swab sterile DNA controlled standard cotton tipped applicator is designed for reliable collection and transport of forensic and DNA evidence. It features a sterile cotton tip and a breathable dry transport tube to help maintain sample integrity.



FS934

SPECIFICATIONS

Construction

- USP grade cotton tip (carded and coiled)
- Polystyrene handle
- Re-closable dry transport tube
- Breathable filter label
- Individually wrapped packaging
- Evidence information label

Materials

- Cotton
- Polystyrene

Performance

- Collects DNA and forensic evidence samples
- Can be used dry or with IPA/sterile water
- Supports sample integrity during transport
- Designed for reaching small or confined areas

Dimensions

- Overall Length: 5.095"
- Tip Diameter: 0.12" tapered to 0.05"
- Tip Length: 0.5"
- Handle Diameter: 0.083"
- Handle Length: 5.032"

Quantity

- 50 per pack

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI