

Swab Protectors

Swab Protectors are designed to cover and protect swabs during transport when immediate air drying is not possible. They help prevent contamination and secure the swab head after collection.



FS936

SPECIFICATIONS

Performance

- Protects swabs from cross-contamination
- Easily slide the protector up the shaft of 6" or 3" swabs

Quantity

- 25 per pack

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI