

Sterile Water 3ml

These 3ml sterile water vials are supplied for controlled dispensing of sterile water during forensic and laboratory procedures. They are commonly used to moisten swabs for improved sample collection.



FT005

SPECIFICATIONS

Construction

- Graduated markings in 0.5ml increments
- Precision dispensing tip
- Easy twist-and-pull top

Capacity

- 3mL

Performance

- Used to moisten swabs for evidence collection
- Enables controlled dispensing by drip or stream

Quantity

- 100 per pack

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI