


Sawtooth “EVIDENCE” Tape

Write-on sawtooth evidence tape is a tamper-indicating tape with a writable strip for sealing and labeling evidence in crime scene and lab settings.



SPECIFICATIONS

Construction

- Dispenser box
- Sawtooth edge design

Quantity

- 5 rolls per case

Performance

- Indicates tampering by tearing and dye smearing
- Easy dispensing with pop-up tape design
- Seals and labels evidence containers

Dimensions

- 1.25" X 108'

Variations

Model No.	Description	Color
FT049	Sawtooth “EVIDENCE” Tape	Blue
FT050	Sawtooth “EVIDENCE” Tape	Clear
FT051	Sawtooth “EVIDENCE” Tape	Green
FT052	Sawtooth “EVIDENCE” Tape	Orange
FT053	Sawtooth “EVIDENCE” Tape	Red
FT054	Sawtooth “EVIDENCE” Tape	Yellow

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI