

C-Arm Accessible Autopsy Carrier

Mopec's C-Arm Accessible Autopsy Carrier is engineered to provide full imaging access during autopsy procedures, allowing seamless use without repositioning. Its open frame design, stability, and smooth mobility support efficient workflow while maintaining excellent visibility and clearance for imaging equipment. Built for precision and durability, it enhances procedural flexibility and operator convenience.



600014

SPECIFICATIONS

Construction

- Meant for use with C-arm type x-ray machines
- (4) locking caster wheels
- Underside of T3603 tray mates to top of the carrier to securely hold the tray in place without the use of mechanical fasteners
- 2" slope for proper drainage, when required

Materials

- 11ga 1.5" SQ 304 Stainless steel frame with a 180-grit polish finish

Dimensions

- 80.5" L X 31" W X 31.8" to 33.9" H (lower to higher end for drainage)

Capacity

- 500 lbs

Compatibility

- T3603 (Plastic Body Tray)

Note:

- Trays sold separately.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI