


# Foot-Elevated Immersion Dissection Table and Dip Tank

The foot-elevated immersion dissection table and dip tank streamlines cadaver handling with durable construction and intuitive operation. Designed for safety and efficiency, it supports the demands of modern anatomy labs and pathology facilities, ensuring reliable performance and ease of use.



## SPECIFICATIONS

### Construction

- Cadaver tray fabricated for strength and easy drainage.
- Stainless steel hinged hood
- 8" casters with brake
- .75" ball valve on drain for hose attachment and fluid control
- Accessory holder (Support sockets) located on all (4) corners to accommodate optional accessories

### Materials

- 304 stainless steel with heliarc welded seams and joints ground and polished to a smooth finish

### Dimensions

- 92.7" L x 30.4" W x 38.7" H

### Capacity

- 300 Lbs

### Performance

- Foot operated hydraulic lifting mechanism on both ends raises immersion tray.
- Automatic lock mechanism locks elevation pedals to keep tray in raised position.

### Options & Accessories

- HO005: Post-Style Leg Holder
- HO006: Post-Style Book Holder
- HO030: Post-Style Instrument Holder

## Related/Alternate/Obsolete Product Model Numbers

Model No.	Description	Status	Notes
600022-01	6" Single Connection	Obsolete	
600022-02	6" Single Connection	Obsolete	
600022-03	6" Single Connection	Obsolete	
600022-04	6" Single Connection	Obsolete	

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI