

Double Scissor Lift Carrier with End/Side Rollers

Mopec’s Double Scissor Lift Carrier with End/Side Rollers is a durable, easy-to-maneuver transport solution designed for safe and discreet body handling. Featuring a double-scissor lift mechanism for smooth height adjustment, it allows loading from either the end or side. The integrated cover provides privacy, while the sturdy construction ensures reliable performance in clinical, forensic, and mortuary settings.



600040-DS-ES-23



600040-DS-ES-27

SPECIFICATIONS

Construction

- End and Side loading
- Fits 23” or 27” trays
- 16” to 56.8” elevation
- Raise or lower the carrier using the handheld controller
- Removable aluminum tubing canopy and black vinyl cover
- (3) heavy duty 6” swivel locking casters
- (1) heavy duty 6” directional locking caster
- (20) Galvanized steel rollers
- Emergency stop button
- (2) locking latch mechanisms at the head and foot of the carrier

Materials

- 11 Gauge square 304 stainless steel carrier frame with a No. 4 polish finish
- Galvanized steel rollers
- HDPE and Nylon omnidirectional rollers
- Removable canopy frame made of aluminum tubing
- Black vinyl removable canopy cover

Dimensions

- **ES-23:** 81.5” L X 26.8” W X 16” to 56.8” H (72.7” H with false frame)
- **ES-27:** 81.5” L X 30.8” W X 16” to 56.8” H (72.7” H with false frame)

Utilities / Power

- Rechargeable battery cord accepts any standard 120V/1hp/60V outlet
- 24V High-capacity battery capable of thrust up to 6000N

Compatibility

- **ES-23:** JC023 and JC023D trays
- **ES-27:** JC027 and JC027D tray
- **Note:** Meant for 3 tier racks, check your application. A 3-tier rack inside a cooler may require a 5-tier lift if the lift is on the floor outside and the rack is inside the cooler

Capacity

- 500 lbs

Note:

- **Trays sold separately.**

Variations

| Model No. | Description | Dimensions |
|-----------------|------------------------------------|--------------------------|
| 600040-DS-ES-23 | Scissor Lift Carrier for 23” Trays | 81.5”L x 26.8”W x 56.8”H |
| 600040-DS-ES-27 | Scissor Lift Carrier for 27” Trays | 81.5”L x 30.8”W x 56.8”H |



www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI