


Hydraulic Covered Cadaver Carrier With Rollers

Mopec carriers are designed to provide safe, secure, and efficient transport of specimens and tissues throughout laboratory, pathology, and mortuary workflows. Built for durability and ease of handling, they help maintain specimen integrity while supporting streamlined operations and compliance with industry standards.



600040-H-27

SPECIFICATIONS

Construction

- Raise or lower the cart using the dual sided mounted foot controls
- Enable fifth wheel steering, neutral, or full brakes with foot controls
- 13" elevation
- (4) 8" diameter wheels
- Removable aluminum tubing canopy and black vinyl cover
- (5) Galvanized steel rollers for smooth transfer of the tray
- (2) locking latch mechanisms at the head and foot of the carrier

Materials

- 11 Gauge square 304 stainless steel carrier frame with a No. 4 polished finish
- Galvanized steel rollers
- The removable canopy frame is made of aluminum tubing
- Hydraulic Stryker base

Dimensions

- 81" L X 31" W X 22" - 35" H (37.4" - 50.4" max height including false frame)

Capacity

- 500 lbs

Compatibility

- JC027 Tray (sold separately)
- JZ029 Cover (included and available separately)

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI