



# Ensemble Portable Vented Trash Can



## SPECIFICATIONS

### Construction

- Provides a vented trash receptacle solution next to your unit's sink
- Includes trash can assembly on casters, flex duct, (2) clamps, and duct stub in side panel of grossing station unit

### Materials

- 304 stainless steel construction with a polished No. 4 finish

### Dimensions

- 15.3" L x 15.2" D x 44.2" H

### Performance

- AX060-L is designed for left-handed Maestro and Maestro Encore grossing stations
- AX060-R is designed for right-handed Maestro and Maestro Encore grossing stations

### Capacity

- 60 lbs

### Compatibility

- Fits all Mopec Maestro and Maestro Encore grossing stations

## NOTE:

- Compatible with all Maestro and Maestro Encore grossing stations.
- Must be added to same side as sink.
- Your lab must have an additional 24" x 24" of space next to your grossing station to accommodate this accessory.

## Variations

Model No.	Description	Dimensions
AX060L	Portable Vented Trash Can, Left	15.3" L x 15.2" D x 44.2" H
AX060R	Portable Vented Trash Can, Right	15.3" L x 15.2" D x 44.2" H

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI