

# 31" Tilting Autopsy Cart

Mopec's autopsy carts are made of stainless steel and are a durable design to support bodies safely and efficiently during examinations. Built with a sturdy frame, options for elevation, and locking casters, it allows staff to transport remains with ease while maintaining a sanitary and controlled work environment.



DB031

## SPECIFICATIONS

### Construction

- For Use with the 31" GA031 tray (sold separately)
- Heavy duty 1.5" stainless steel tubing frame
- Heavy duty 8" casters with total-lock feature
- Single-sided hand crank with 7-1/2" lift adjustment

### Materials

- 1.5" Square 304 stainless steel construction with a 180-grit finish

### Dimensions

- 68.1" L x 31.3" W x 36"-42.5" H (lower to higher end when elevated)

### Capacity

- 500 lbs

### Compatibility

- GA031 tray with sloped sides (sold separately)
- Compatible with Mopec's CA CB, CC autopsy sinks, but does not work with CC200 due to elevation

### Note:

- Trays sold separately.

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI