

Cartridge Tube Refills, Various Colors

These Cartridge Tube Refills are replacement cartridge tubes of forensic casting material designed for use with the Automix Dispensing System. Made from a non-toxic polyvinylsiloxane evidence casting material that is flexible, highly resistant to tearing, and provides strong dimensional stability for accurate evidence casting.



For reference only: Casting Starter Kit

SPECIFICATIONS

Construction

- Non-toxic polyvinylsiloxane evidence casting material provided in various colors
- To be loaded into the Dispensing Gun

Options & Accessories

- FS371: Casting Starter Kit
- FS372: Dispensing Gun

Quantity

- 4 per package

Capacity

- 75 mL

Performance

- Cures in 4 minutes at 68°F
- Brown: reduce glares, ideal for use in casting gun barrels and shell casings
- Transparent: no distortion, ideal for use on rough or curved surfaces
- White: lifts and preserve latent prints from textured surfaces, ideal for use in casting gun barrels and palm prints

Note:

- Photo for reference only. Cartridge Tube Refills included only.

Variations

Model No.	Color
FS375	Brown Cartridge Tube Refills
FS376	Transparent Cartridge Tube Refills
FS377	White Cartridge Tube Refills

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI