

Casting Putty Kit, Various Colors

This Casting Putty Kit is used to create highly detailed casts of tool marks and to lift latent prints from rough or curved surfaces. It comes in various colors, mixes quickly, sets in minutes, and can produce precise impressions even in cold temperatures.



SPECIFICATIONS

Construction

- Kit consists of:
- (1) Mixing Card
- (1) Spatula
- (1) Tube of Catalyst
- (1) Casting Material
- (1) Tube of Casting Putty
- Used when highly detailed castings of a tool and jimmy marks are needed or for lifting dusted latent prints off rough or round surfaces
- Excellent at rendering small details
- Short curing time
- High contrast for microscopic observations
- Works in below freezing temperatures

Options and Accessories

- FS396: Additional Tube of Catalyst

Performance

- Mix the putty and catalyst on a mixing card with the spatula and apply it to the surface being cast
- Cast should set between 6-8 minutes
- In below freezing temperatures, cast should be set within 15 minutes
- Black:
- Use to lift dusted latent fingerprints from rough or textured surfaces where regular lift tape will not work
- Use black casting putty with white, silver, or fluorescent

powders

- White:
- Use to lift dusted latent fingerprints from rough or textured surfaces where regular lift tape will not work
- Use white casting putty with black powders
- Gray:
- Provides contrast between light-colored surfaces and shadows produced from side-lighting
- Brown:
- Dark color increases the visibility of detail in the casts

Variations

Model No.	Description	Color
FS392	Casting Putty Kit	Black
FS393	Casting Putty Kit	White
FS394	Casting Putty Kit	Gray
FS395	Casting Putty Kit	Brown

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI