

Chlorate/Bromate Test Kits

The Chlorate / Bromate Test Kit provides fast, safe, and reliable on-the-spot detection of chlorate ions using sealed, contact-free ampoule technology. It delivers clear results in seconds, requires no mixing or measuring, and features a built-in neutralization step for safe handling.



FS410

SPECIFICATIONS

Construction

- Easy-to-crush numbered, glass ampoules

Performance

- 3-year expiration date
- Specific and stable test of Chlorate ion (an oxidizing ingredient of low explosives)
- Self-contained neutralization reagent of strong acid is used after the test

- Clearly visible colors for identification
- Immediate, reliable results
- Environmentally safe
- No chemical leakage
- Chemical-free experience for operator
- High sensitivity to trace amounts
- Simple to use

Materials

- Plastic
- Glass

Quantity

- 10 test kits

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI