

CS Forensic Photo Light System

The CS Forensic Photo Light System is a rechargeable, waterproof light with multiple UV and blue wavelengths for detecting biological fluids and trace evidence. It comes as a complete kit with orange and yellow goggles, filters, charger, and carrying accessories, making it ready for fast, reliable forensic work in the field or lab.



FS449

SPECIFICATIONS

Construction

- Includes:
- (1) CS Scan Light
- (1) Orange Goggle
- (1) Yellow Goggle
- (1) Charger
- (1) Adapter
- (1) Carry Strap
- (1) Case
- (1) CS Camera Filter System
- Rectangular anti-roll head
- Beam angle of 8 degrees
- Rechargeable Lithium-Ion Battery System
- Offers various color and wavelengths in 1 light:

- LED colors
UV, blue, and white
- Wavelength options
365, 395, 450, and 470nm
White light is also available

Dimensions

- Head: 3.8" diameter
- Handle:
• 2.1" diameter
• 11.1" Length

Materials

- Aluminum 6061-T6 hard anodized housing
- TPR bumpers
- Glass lens

Performance

- Helps speed up forensic searches and reduces cross-contamination at the scene
- Trusted to look for body fluids, accelerants, bite marks, bruising, and trace evidence
- Waterproof
- Fire resistant, meeting the requirements of NFPA 1971-8.6 (2013)
- CE-EMC certification
- 4 hour recharge time
- 6 hour battery life
- Can be easily cleaned and disinfected after each use

Note:

- **Limited Lifetime Warranty**

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI