

# UV Forensic Light System

The UV Forensic Light System is a compact, high-output illumination tool designed to enhance the detection of biological fluids, trace evidence, and latent residues for convenient field or lab use. The light source has a fixed wavelength output of 365nm and features a rechargeable design for extended use without downtime.



FS450

## SPECIFICATIONS

### Construction

- Forensic light source  
365nm  
Rectangular anti-roll head  
Rechargeable Lithium-Ion Battery System
- Also includes:
  - (1) Yellow goggles
  - (1) Charger
  - (1) Adapter
  - (1) Carry strap
  - (1) Case

### Performance

- Trusted to look for body fluids, accelerants, bite marks, bruising, and trace evidence
- Forensic evidence is visible up to 20ft away
- High quality and long lasting
- Withstands physical impact and is waterproof
- Accurate to +5m which helps visualize fluorescence and absorption in evidence
- Can be easily cleaned and disinfected after each use

### Materials

- Aluminum

### Options & Accessories

- FS451: Blue Rechargeable Forensic Light Source, 365nm
- FS452: Cyan Rechargeable Forensic Light Source, 495nm
- FS453: Green Rechargeable Forensic Light Source, 525nm
- FS454: All-In-One, 470, 495, 525, and 850nm

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI