

Hot Plate, 120v

Hot plate with a 3" heating surface is used with aluminum dishes inside a fuming chamber to accelerate cyanoacrylate fuming. It features an easy on/off switch and is compatible with 2" diameter aluminum trays for latent print development.



FS502

SPECIFICATIONS

Construction

- Easy on/off switch

Option and Accessories

- FS501: Aluminum Tray

Performance

- Accelerate the cyanoacrylate fuming process by using an aluminum tray (FS501, sold separately) and aluminum dish inside a fuming chamber

Materials

- Aluminum

Dimensions

- 2" diameter X 1.5" deep

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI