

Ninhydrin Print Fixative

The Ninhydrin Print Fixative is used to permanently preserve prints developed using Ninhydrin. Simple application of spraying the fixative over the print and letting it air dry.



FS526

SPECIFICATIONS

Construction

- Permanently lock in and preserve prints developed with Ninhydrin

Performance

- Spray over the print and let it air-dry

Capacity

- 8oz

Materials

- Ninhydrin

Options & Accessories

- FS527: Ninhydrin Print Remover

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI