


Back-Up Hinged Latent Print Lifters, Various Sizes

A hinged latent print lifter designed to improve collection and documentation by incorporating a pre-attached backing card with a preprinted reverse side for case information, while allowing the lifter to be opened for backlight (transmitted light) examination of lifted prints.



SPECIFICATIONS

Construction

- Standard adhesive-backed hinged lifter design
- Preprinted reverse side for case documentation
- White or black high-gloss backing card affixed to lifter
- .012" thick backing card material

Performance

- Lifts and preserves latent fingerprints and palm prints
- Provides immediate background contrast for print visibility
- Enables backlit (transmitted light) examination in lab
- Improves detail visibility during photography and analysis
- Ensures case information is recorded on backing card

Quantity

- 24 per pack: FS598, FS599, FS602, and FS603
- 12 per pack: FS600, FS601, FS604, FS605

Variations

Model No.	Description	Dimensions
FS598	Black Back-Up Hinged Latent Print Lifters	1.5" X 2"
FS599	White Back-Up Hinged Latent Print Lifters	1.5" X 2"
FS600	Black Back-Up Hinged Latent Print Lifters	2" X 4"
FS601	White Back-Up Hinged Latent Print Lifters	2" X 4"
FS602	Black Back-Up Hinged Latent Print Lifters	4" X 4"
FS603	White Back-Up Hinged Latent Print Lifters	4" X 4"
FS604	Black Back-Up Hinged Latent Print Lifters	6" X 9"
FS605	White Back-Up Hinged Latent Print Lifters	6" X 9"

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI