

Luminol Reagent

A chemiluminescent blood detection reagent used for locating trace amounts of blood over large search areas.



SPECIFICATIONS

Construction

- Field kits include spray bottle, liquid reagent, pre-weighed Luminol powder, and instructions

Performance

- Detects trace blood evidence
- Suitable for large-area crime scene searches
- Uses chemiluminescent Luminol method
- Field-applicable blood detection system

Volume

- FS755: 2oz
- FS756: 8oz
- FS757: 16oz
- FS758: 4oz X 6
- FS759: 5g

Variations

Model No.	Description
FS755	Luminol Field Kit
FS756	Luminol Field Kit
FS757	Luminol Field Kit
FS758	Luminol Field Kit, Powder Packs
FS759	Luminol Powder

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI