

Luminol Reagent

The Luminol Reagent is formulated for use in the chemiluminescent detection of latent blood evidence. When applied, it reacts with trace amounts of blood to produce a visible glow, allowing investigators to locate otherwise invisible stains.



SPECIFICATIONS

Construction

- Luminol powder is available for separate purchase
- Kits are available in various sizes
- Kits include:
 - (1) Spray Bottle
 - (1) Measured amount of liquid reagent
 - (1) Pre-weighed vial of Luminol powder
 - (1) Complete Instructions

Performance

- Detects trace blood suitable for large-area crime scene searches
- Presumptive test for blood
 - Does not discriminate between human and non-human species of blood
- Uses the chemiluminescent Luminol method

Capacity

- FS755: 2oz
- FS756: 8oz
- FS757: 16oz
- FS578: (6) 4oz packs
- FS759: 5g

Variations

Model No.	Description	Capacity
FS755	Luminol Field Kit	2oz
FS756	Luminol Field Kit	8oz
FS757	Luminol Field Kit	16oz
FS758	Luminol Field Kit, Powder Packs	(6) 4oz packs
FS759	Luminol Powder	5g

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI