

Ardrox Forensic Reagent

A fluorescent liquid dye used to enhance latent fingerprints for examination with lasers, forensic light sources, and UV lamps.



SPECIFICATIONS

Performance

- Enhances latent fingerprints
- Fluoresces under UV and forensic light sources
- Optimized viewing between 450–480 nm
- Can be used full-strength or diluted
- One-year shelf life

Variations

Model No.	Description	Volume
FS760	Ardrox Concentrate	6 X 20mL
FS761	Ardrox Forensic Reagent	500mL

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI