

Ardrox Forensic Reagent

A fluorescent liquid dye used to develop prints for examination with lasers, forensic light sources, and UV lamps. Developed prints fluoresce under 365nm light, but better results are found when viewing between 450–480nm.



SPECIFICATIONS

Construction

- Fluorescent liquid dye
- Used to develop prints
- Can be used full-strength or diluted with ethanol or methanol
- Shelf-life of 1 year

Performance

- Fluoresces under lasers, forensic light sources, and UV lamps
- Prints fluoresce under lights in 365nm, but better results are found when viewing between 450–480nm

Materials

- Ardrex

Capacity

- FS760: (6) 20mL
- FS761: 500mL

Variations

Model No.	Description	Capacity
FS760	Ardrex Concentrate	(6) 20mL
FS761	Ardrex Forensic Reagent	500mL

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI