

# Basic Yellow

A fluorescent yellow dye used after cyanoacrylate development to enhance latent fingerprints for examination with forensic light sources or UV lamps.



## SPECIFICATIONS

### Performance

- Enhances cyanoacrylate-developed latent prints
- Fluoresces between 365–485 nm
- Used with forensic light sources or UV lamps
- Improves latent print visibility and contrast

### Variations

Model No.	Description	Capacity
FS769	Basic Yellow Powder	25g
FS770	Basic Yellow Reagent	1L

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI