

Leuco Crystal Violet

A chemical reagent used to identify and enhance bloodstains and bloody latent prints by reacting with the heme component of blood to produce a violet coloration.



SPECIFICATIONS

Performance

- Enhances bloody latent prints
- Identifies non-visible bloodstains
- Increases contrast of visible blood patterns
- Reacts with heme in blood
- Commonly used on porous surfaces

Construction

- Available as powder or field kit formulations

Variations

| Model No. | Description | Capacity |
|-----------|--------------------------------|----------|
| FS775 | Leuco Crystal Violet Field Kit | 16 oz |
| FS776 | Leuco Crystal Violet Field Kit | 32 oz |
| FS777 | Leuco Crystal Violet Powder | 25 g |

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI