

# RAY Reagent

A fluorescent dye mixture used to enhance cyanoacrylate-developed latent fingerprints, especially on plastic and non-porous surfaces.



## SPECIFICATIONS

### Performance

- Enhances cyanoacrylate (CA)-developed prints
- Particularly effective on plastic bags and similar surfaces
- Fluoresces under forensic light sources
- Visualized between 415–530 nm
- Improves contrast and ridge detail

### Capacity

- 1 L

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI