

# Rhodamine 6G Solution

A fluorescent dye solution used to enhance cyanoacrylate-developed latent fingerprints for examination with lasers or forensic light sources.



## SPECIFICATIONS

### Performance

- Enhances cyanoacrylate (CA)-developed latent prints
- Improves visibility on non-porous surfaces
- Fluoresces under laser or forensic light source
- Provides strong ridge detail contrast

### Variations

Model No.	Description	Capacity
FS782	1,2-Indandione Powder	1L
FS783	1,2-Indandione Powder	500mL

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI