

Postmortem Specimen Collection Kit

The Postmortem Specimen Collection Kit is designed for safely collecting and transporting biological specimens to medical examiners or crime laboratories for analysis. It meets all shipping regulations and includes pre-labeled, secure components for reliable and compliant specimen handling.



SPECIFICATIONS

Construction

- (1) Tamper-evident, pre-sealed Kit Box with biohazard label affixed
- (1) Kit Instruction Sheet/FDA Insert
- (1) Specimen Bottle Holder
- (1) 20 ml Peripheral Blood Vial containing 400 mg of sodium fluoride powder with 20 ml "Fill To Mark"
- (1) 20 ml Aorta Blood Vial containing 400 mg of sodium fluoride powder
- (1) 20 ml "Other" Vial
- (1) 90 ml Liver Specimen Bottle with label affixed

- (1) Resealable Bag, containing a liquid-absorbing sheet
- (1) Kit Box Shipping Seal for resealing kit box after specimen collection

Performance

- Collects and secures biological specimens for analysis
- Ensures compliance with shipping and transport regulations
- Provides organized, pre-labeled components for accuracy
- Maintains specimen integrity during transport

Dimensions

- Box: 6"L x 4.75"W x 2.5"H
- Bag: 8" x 8"

Quantity

- 25 kits per case

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

800.362.8491
800 Tech Row, Madison Heights, MI