

Postmortem DWI Specimen Collection Kit

The Postmortem DWI Specimen Collection Kit is designed for the secure collection of blood and urine samples from deceased individuals for toxicology analysis. It includes compliant packaging and tamper-evident components to ensure safe and regulated transport to laboratories or medical examiners.



FS804

SPECIFICATIONS

Construction

- (1) Tamper-evident Pre-sealed Kit Box with biohazard label affixed
- (1) FDA Product Insert
- (1) Inner Box Blood Tube and Urine Bottle Holder with security flap
- (1) Tamper-evident Evidence Seal for sealing inner kit box
- (1) 100 ml Urine Specimen Bottle
- (2) 10 ml Blood Collection Tubes (gray-stoppered) containing 100 mg of sodium fluoride and 20 mg of potassium oxalate
- (1) 16 g x 4" Needle
- (1) Resealable Bag containing a liquid-absorbing sheet

- (3) Tamper-evident Specimen Security/Identification Seals
- (1) Tamper-evident Kit Box Shipping Seal for resealing outer kit box after specimen collection

Performance

- Collects blood and urine samples for toxicology testing
- Ensures secure, tamper-evident specimen handling
- Meets shipping regulations for biological materials
- Maintains sample integrity during transport

Dimensions

- Box: 6" L x 4.75" W x 2.5"H

Quantity

- 25 kits per case

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

800.362.8491
800 Tech Row, Madison Heights, MI