

Precut Latent Print Lifting Tape Pad

This precut transparent lifting tape comes in a convenient peel-off pad for quick and easy fingerprint collection. Each strip is ready to use, allowing for efficient lifting and transfer to backing cards.



FS898 and FS899

SPECIFICATIONS

Construction

- Peel-off pad
- Precut strips
- White leader tabs

Performance

- Lifts developed fingerprints cleanly and clearly
- Peel-off design allows quick, efficient use
- Transfers prints securely to backing cards

Quantity

- 25 strips per pad

Variations

| Model No. | Dimensions |
|-----------|------------|
| FS898 | 2" x 6" |
| FS899 | 4" x 6" |

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI