

Latent Print Lifting Tape

The Latent Print Lifting Tape is designed to lift prints from a variety of sources, featuring uniform adhesive application for consistent results. Available in multiple widths, finishes, and frosted or transparent options. They are packaged to prevent contamination in resealable bags or plastic tubes.



SPECIFICATIONS

Construction

- Available in resealable bags or plastic tubes to prevent contamination
- Available in frosted and transparent finishes
- 1" cores

Performance

- Uniform application of adhesive and are heat treated and pressure wound to eliminate imperfections like fish-eyes
- Lifts fingerprints and trace evidence cleanly

Dimensions

- All variations have a length of 360"
- Width sizes include: 1.5", 2", 3", and 4"

Variations

Model No.	Description	Dimensions
FS900	Frosted in Bag Latent Print Lifting Tape	360" L x 1.5" W
FS901	Frosted in Bag Latent Print Lifting Tape	360" L x 2" W
FS902	Frosted in Bag Latent Print Lifting Tape	360" L x 3" W
FS903	Frosted in Bag Latent Print Lifting Tape	360" L x 4" W
FS904	Frosted in Tube Latent Print Lifting Tape	360" L x 1.5" W
FS905	Frosted in Tube Latent Print Lifting Tape	360" L x 2" W
FS906	Frosted in Tube Latent Print Lifting Tape	360" L x 3" W
FS907	Frosted in Tube Latent Print Lifting Tape	360" L x 4" W
FS908	Transparent in Bag Latent Print Lifting Tape	360" L x 1.5" W
FS909	Transparent in Bag Latent Print Lifting Tape	360" L x 2" W
FS910	Transparent in Bag Latent Print Lifting Tape	360" L x 3" W
FS911	Transparent in Bag Latent Print Lifting Tape	360" L x 4" W
FS912	Transparent in Tube Latent Print Lifting Tape	360" L x 1.5" W
FS913	Transparent in Tube Latent Print Lifting Tape	360" L x 2" W
FS914	Transparent in Tube Latent Print Lifting Tape	360" L x 3" W
FS915	Transparent in Tube Latent Print Lifting Tape	360" L x 4" W

MOPEC
GROUP



www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI