

Inking Slabs

These Inking Slabs provide a smooth, durable surface for controlled fingerprint ink application and are available in polished stainless steel or 0.25" thick glass. They are designed for stable use in fingerprinting stations and laboratory settings.



SPECIFICATIONS

Construction

- Stainless steel slab
Pre-drilled for tabletop mounting
- Glass slab
0.25" thick
Provided with non-skid rubber feet

Options and Accessories

- 0.25" thick wood backing plate to raise the slab off of the table for easy cleaning (FT101 and FT102)

Performance

- Provides controlled surface for applying fingerprint ink
- Available in multiple sizes for different applications

Materials

- Glass
- Rubber
- Wood
- Stainless steel

Variations

Model No.	Description	Dimensions
FT097	Glass Inking Slab	4" X 10"
FT098	Glass Inking Slab	5" X 12"
FT099	Stainless Steel Inking Slab	4" X 10"
FT100	Stainless Steel Inking Slab	5" X 12"
FT101	Wood Riser/Backing Plate	4" X 10"
FT102	Wood Riser/Backing Plate	5" X 12"

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI