

Presumptive Drug Tests Portable Kit

This portable presumptive drug identification kit is designed for field narcotics operations where mobility and ease of use are essential. It presumptively identifies up to 25 controlled substances for rapid on-site screening.



FT149

SPECIFICATIONS

Construction

- Includes:
 - (2 ea) Tests A, B, C, D, E, G, J, K, L, and U
 - (1) Desk Chart
 - (1) Polytesting Instruction Manual
 - (1) Polytesting Report Pad
 - (1) Plastic Carrying Case (10.5" W x 9" H x 3" D)

Performance

- Portable presumptive drug screening
- Identifies up to 25 controlled substances
- Designed for field narcotics operations
- Rapid on-site testing
- Long shelf-life reagents (proper storage)

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI