

# Rapid Drug Detection Dip Sticks

A simple, on-site urine screening tool designed to quickly detect the presence of synthetic cannabinoid compounds (such as K2/Spice), using a dip-and-read method that delivers results in about five minutes; it provides qualitative detection at a specific cutoff level and is intended for fast, efficient screening in forensic or controlled testing environment



FT263

## SPECIFICATIONS

### Construction

- Simple and easy procedure: Collect, dip, and read results
- Tests urine for the presence of synthetic cannabinoid compounds Found to be several times stronger than marijuana
- Qualitative detection of JWH-018 and JWH-073 metabolites in urine at 30 ng/mL cutoff level

### Quantity

- 25 per pack

### Performance

- For on-site testing
- On-site test results in 5 min.
- Average window of detection is 72 hours following a single low dose; in case of chronic use the window is much longer

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI