

Polybag Evidence Tubing Dispenser, Various Sizes

A sturdy rack-style tool used in forensic and laboratory packaging to hold rolls of polybag tubing and dispense them in a controlled way, making it easier to pull and cut custom-length bags for evidence items that don't fit standard containers; it helps keep the tubing organized, prevents tangling or unraveling, and improves efficiency when preparing evidence packaging in controlled environments.



SPECIFICATIONS

Construction

- Designed to dispense polybag evidence tubes
- Make polybags easily accessible
- Available in 2 sizes

Options & Accessories

- Polybag evidence tubing, various sizes
 FT291: 3" x 1125'
 FT292: 7" x 1125'
 FT293: 12" x 900'
 FT294: 24" x 645'

Performance

- Dispense any length of tubing needed to package evidence

Note:

- Polybag evidence tubing is sold separately

Variations

Model No.	Description	Dimensions
FT289	Polybag Evidence Tubing Dispenser	12"
FT290	Polybag Evidence Tubing Dispenser	24"

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI