

Polybag Evidence Tubing, Various Sizes

A forensic packaging device designed to heat-seal lengths of polybag tubing used for securing and storing evidence, allowing users to create custom-sized, tamper-evident containment for items collected at a scene; it includes a built-in cutter for trimming the tubing cleanly and is intended to streamline the sealing process for crime scene and laboratory evidence handling.



SPECIFICATIONS

Construction

- Evidence collection bags
- Available in 4 sizes

Materials

- 4 mil

Options & Accessories

- Polybag evidence tube dispenser, various sizes
FT289: 12”
FT290: 24”
- Polybag evidence tube sealer, various sizes
FT295: 14”
FT296: 20”
FT297: 26”

Performance

- Package evidence that will not fit in commonly-sized bags

Note:

- Polybag Evidence Tubing Dispenser is sold separately

Variations

Model No.	Description	Dimensions
FT291	Polybag Evidence Tubing	3" x 1125'
FT292	Polybag Evidence Tubing	7" x 1125'
FT293	Polybag Evidence Tubing	12" x 900'
FT294	Polybag Evidence Tubing	24" x 645'

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI