

Polybag Evidence Tube Sealer, Various Sizes

A forensic packaging device designed to heat-seal lengths of polybag tubing used for securing and storing evidence, allowing users to create custom-sized, tamper-evident containment for items collected at a scene; it includes a built-in cutter for trimming the tubing cleanly and is intended to streamline the sealing process for crime scene and laboratory evidence handling.



SPECIFICATIONS

Construction

- Designed to dispense polybag evidence tubes and seal evidence
- Set of rollers
- Cutter

Options & Accessories

- Polybag evidence tubing, various sizes
 - FT291: 3" x 1125'
 - FT292: 7" x 1125'
 - FT293: 12" x 900'
 - FT294: 24" x 645'

Performance

- Can seal any length up to 26" W
 - Roll of polybag rides on a set of rollers to dispense any length of tubing needed to securely enclose the evidence
- Heat-seal one end of the tubing, turn it around and seal the other end

Note:

- Polybag evidence tubing is sold separately

Variations

Model No.	Description	Dimensions
FT295	Polybag Evidence Tube Sealer w/Cutter	14"
FT296	Polybag Evidence Tube Sealer w/Cutter	20"
FT297	Polybag Evidence Tube Sealer w/Cutter	26"

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI