

Rubber/Gelatin Footprint Lifters

Forensic tools designed to recover and preserve developed fingerprints from a variety of textured or irregular surfaces—such as paper, cardboard, or wallpaper—without damaging the evidence; they use a low-adhesive gelatin layer on a flexible backing with a hinged cover to securely lift and protect prints for examination and storage.



SPECIFICATIONS

Construction

- Designed specifically for lifting footprints
- Different colors available

Quantity

- 2 per pack

Materials

- Rubber backing
- Thick adhesive gelatin

Dimensions

- 5" x 14"

Performance

- Easy application and removal to any surface
- Useful for collecting minute particles of dust or dirt left by footwear on rough surfaces
- Lifted footprints fade over time
Lower the temperature, the slower the fading
- Oblique lighting will show details not visible before photography

Variations

Model No.	Description	Colors
FT301	Rubber/Gelatin Footprint Lifter	Black
FT302	Rubber/Gelatin Footprint Lifter	White
FT303	Rubber/Gelatin Footprint Lifter	Transparent

MOPEC
GROUP



www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI