

Rubber/Gelatin Hinged Latent Print Lifters

Forensic tools designed to recover and preserve developed fingerprints from a variety of textured or irregular surfaces—such as paper, cardboard, or wallpaper—without damaging the evidence; they use a low-adhesive gelatin layer on a flexible backing with a hinged cover to securely lift and protect prints for examination and storage.



SPECIFICATIONS

Construction

- Designed specifically for lifting fingerprints
- Different colors and sizes available

Quantity

- 12 per pack

Materials

- Rubber backing
- Low-adhesive gelatin
- Clear mylar cover sheet

Performance

- Easy application and removal to any surface
- Can collect developed prints from rough, grained, or irregular surfaces

Variations

Model No.	Description	Dimensions	Colors
FT301	Rubber/Gelatin Hinged Latent Print Lifters	2" X 2"	Black
FT302	Rubber/Gelatin Hinged Latent Print Lifters	2" X 2"	White
FT303	Rubber/Gelatin Hinged Latent Print Lifters	2" X 2"	Transparent
FT304	Rubber/Gelatin Hinged Latent Print Lifters	2" X 4"	Black
FT305	Rubber/Gelatin Hinged Latent Print Lifters	2" X 4"	White
FT306	Rubber/Gelatin Hinged Latent Print Lifters	2" X 4"	Transparent
FT307	Rubber/Gelatin Hinged Latent Print Lifters	4" X 4"	Black
FT308	Rubber/Gelatin Hinged Latent Print Lifters	4" X 4"	White
FT309	Rubber/Gelatin Hinged Latent Print Lifters	4" X 4"	Transparent

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI