

Blood and Urine Specimen Collection Kit

A kit with a comprehensive forensic solution designed to ensure that both blood and urine samples are properly collected, preserved, and clearly identified for laboratory analysis, helping maintain accuracy and chain-of-custody throughout the process; it also includes detailed instructions to guide users and reduce errors during collection.



FS310

SPECIFICATIONS

Construction

- Kit is comprised of:
- (1) Pre-sealed Kit Box with bio-hazard label affixed
- (1) Kit Instruction Sheet / FDA Insert
- (1) Subject Consent Form
- (1) Inner Box Blood Tube Holder and Urine Bottle Holder with Investigating Officer's Report and Chain of Possession printed on security flap
- (1) Tamper-evident Evidence Seal for sealing inner box
- (2) 10 ml Blood Collection Tubes (gray-stoppered) containing 100 mg of sodium fluoride and 20 mg of potassium oxalate

- (1) Multi-sample Needle
- (1) Safety Tube and Needle Holder
- (1) PVP Prep Pad (alcohol-free)
- (1) Resealable Bag containing a liquid-absorbing sheet
- (3) Tamper-evident Specimen Security / Identification Seals
- (1) Tamper-evident Kit Box Shipping Seal for resealing outer kit box after specimen collection

Dimensions

- Pre-sealed kit box
6" L x 4.75" W x 2.5" H
- Multi-sample Needle
20 g x 1.5"

Quantity

- 25 kits

Performance

- Ensure that crime laboratory personnel will receive secure, undamaged, properly collected, and clearly labeled blood specimens for analysis
- Step-by-step instructions included to streamline collection and clarify responsibilities
Defines and separates the responsibilities of the investigating officer and medical personnel
- Boxed kit complies with U.S. Postal and commercial shipping regulations for biological specimens if mailed

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI