

DNA Blood Stain Specimen Collection Kit

The DNA Bloodstain Collection Kit is available in envelope or box configurations to securely collect samples and protect the chain of custody at crime scenes. Each case includes 25 kits pre-packed with sterile swabs, clear instructions, and tamper-evident seals to ensure contamination-free specimens for lab analysis.



FT312



FT313

SPECIFICATIONS

Construction

- Kit is comprised of:
- FT312: (1) Tamper-evident, pre-sealed kit envelope
- FT313: (1) Kit shipping box
- (1) Kit instruction sheet / FDA insert
- Pair of nitrile barrier gloves
- Distilled water dispenser for collection of dry blood samples

- Police evidence seal for resealing kit envelope after specimen collection
- (5) Sterile, cotton-tipped 3" swabs
- (5) Swab boxes
- (5) Swab box envelopes to prevent cross-contamination

- 25 kits

Quantity

- Sterile, cotton-tipped 3" swabs 2 per package

Variations

Model No.	Description	Version	Dimensions
FT312	DNA Blood Stain Specimen Collection Kit	Envelope	6" L x 9" W
FT313	DNA Blood Stain Specimen Collection Kit	Box	6" L x 4.75" W x 2.5" H

MOPEC
GROUP



www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI