

Suspect DNA Fingertick Collection Kit

A DNA collection solution kit designed for obtaining blood samples in forensic investigations through a simple finger-prick method, allowing use by non-medical personnel in many jurisdictions; the collected sample is preserved on a card, then sealed and sent to a crime laboratory for DNA analysis and identification.



FT314

SPECIFICATIONS

Construction

- Kit is comprised of:
- (1) Tamper-evident, pre-sealed kit envelope
- (1) Kit instruction sheet / FDA insert
- (1) Fingertick blood information card containing a blood collection card
- (1) Pair of nitrile barrier gloves
- (1) Alcohol prep pad
- (1) Sterile gauze pad
- (1) 21 gauge fingerstick lancet
- (1) Adhesive bandage

- (1) Resealable bag containing a desiccant packet
- (1) Tamper-evident police evidence seal
- (1) Kit mailing envelope
- (1) Kit shipping seal

Dimensions

- Tamper-evident, pre-sealed kit envelope
6" L x 9" W
- Kit mailing envelope
7.5" L x 10.5" W

Quantity

- 25 kits

Performance

- Designed to be used by non-medical personnel
- How to use:
Use the sterile lancet provided
Prick the subject's left or right ring finger
Squeeze the finger allowing the blood to fill the inner circle on the blood collection card
- Kit is sealed and then forwarded to the crime laboratory for analysis

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI