

Suspect DNA Buccal Swab Collection Kit

A non-invasive DNA collection solution kit designed for use by non-medical personnel, allowing cheek cell samples to be easily obtained and then securely transported to a crime laboratory for identification and comparison with evidence samples; it is part of a suspect DNA collection system



FT315

SPECIFICATIONS

Construction

- Kit is comprised of:
- (1) Tamper-evident, pre-sealed kit envelope
- (1) Kit instruction sheet / FDA insert
- (1) Buccal swab information card containing a swab envelope
- (1) Resealable bag containing a desiccant packet
- (1) Pair of nitrile barrier gloves
- (1) Sterile, cotton-tipped swabs
- (2) Tamper-evident police evidence seals
- (1) Kit mailing envelope
- (1) Kit shipping seal

Dimensions

- Tamper-evident, pre-sealed kit envelope
6" L x 9" W
- Kit mailing envelope
7.5" L x 10.5" W

Quantity

- Sterile, cotton-tipped swabs
2 per package
- 25 kits

Performance

- Least intrusive out of all the other suspect DNA specimen collection kits
- Designed to be used by non-medical personnel
- Cotton-tipped swabs are placed in the subject's mouth, the cheek is thoroughly swabbed, resulting in the transfer of buccal cells onto the cotton bulb of the swab
- Kit can be either mailed or hand-delivered to the crime laboratory for "known" DNA analysis and comparison with "unknown" DNA from a crime scene

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI