

Urine Specimen Collection Kits

These kits are a pre-packaged solution designed for secure and reliable urine sample collection in forensic and DWI testing situations, helping maintain proper handling, identification, and chain-of-custody procedures; it is organized for efficiency and consistency in the field.



FT324



FT325

SPECIFICATIONS

Construction

- Kit is comprised of:
- (1) Pre-sealed kit box with biohazard label affixed
- (1) Kit instruction sheet / FDA insert
- (1) Specimen bottle holder
- (1-2) 100 ml Urine specimen bottle(s)
- (1) Resealable bag containing a liquid-absorbing sheet

- (1) Investigating officer's report and chain of possession label for the resealable bag
- (1-2) Tamper-evident specimen security / identification seal(s)
- (1) Tamper-evident kit box shipping seal for resealing kit box after specimen collection

Quantity

- 25 kits

Dimensions

- Pre-sealed kit box:
5" L x 3.75" W x 2.5" H

Materials

- Specimen bottle holder
- Polyurethane foam

Performance

- Specifically designed for driving while intoxicated (DWI) specimen collection

Variations

Model No.	Description	Number of Bottles	Number of Seals
FT324	Urine Specimen Collection Kit	1	1
FT325	Urine Specimen Collection Kit	2	2

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI