

4.5X Fingerprint Classifier Magnifier

A professional-grade optical tool designed for detailed fingerprint analysis and classification. It uses a 4.5× flat-field lens that provides sharp, distortion-free viewing across the entire image, making it easier to examine fine ridge details like bifurcations and ridge endings.



FT328

SPECIFICATIONS

Construction

- 4.5X magnifier lens
- Finest optics, rated highest in:
 - Center sharpness
 - Overall sharpness
 - Brightness
 - Contrast
 - Color correction
 - Magnification
 - Ease of handling

Materials

- Flat-field lens

Performance

- Distortion-free
- View is crystal clear and sharply focused across the entire field
- For fingerprint analysis
- An independent study by the Science & Practices Committee of the California State Division of the International Association of Identification rated this magnifier the best overall of the seven top models tested

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI