

4X Dome Magnifiers

A handheld optical tool designed for detailed examination of forensic evidence, such as fingerprint record cards and fluorescent fingerprint powders. With 10× magnification, it provides a clear, enlarged view of fine details critical for analysis.



SPECIFICATIONS

Construction

- 4X magnifier lens
- Dome-shaped
- Lightweight and unbreakable

Materials

- Lint-free, velvet-like drawstring bag (included)
- Solid acrylic

Performance

- Scratch-resistant
- Provide distortion-free, crystal clear magnification across the entire magnifier field

Variations

Model No.	Description	Diameter
FT329	4X Dome Magnifier	2"
FT330	4X Dome Magnifier	2.5"
FT331	4X Dome Magnifier	3"

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI